

Gewinde-Adapter 3/4" LaboPlast^{®}

Für StopCock-Ablasshähne und Kugelhähne.

Gewinde auüen	Gewinde 2 auüen	Farbe	Material	Menge pro VE	Bestell Nr.
22x4 (Hünnersdorff)	G 3/4"	rot	PP	1	6.089 262
20x4 (Kautex)	G 3/4"	blau	PP	1	9.011 095
G 3/4"	G 3/4"	schwarz	PP	1	6.078 922
R 3/4"	G 3/4"	grün	PP	1	9.011 096
S56x4 (Tri-Sure)	G 3/4"	grün/weiß	PP/PE	1	7.608 044
S70x6 (Mauser)	G 3/4"	grün/weiß	PP/PE	1	7.608 086
G 2"	G 3/4"	grün/weiß	PP/PE	1	6.901 001