

Splash head adapter, borosilicate glass 3.3, straight

Suitable to be used in distillation assemblies to prevent raw liquid step over from flask to condenser.

Material: Boro 3.3

Cone[s] straight	Socket[s] straight	PK	Cat. No.
NS29/32	NS29/32	1	6.255 503

Zubehör:

PTFE sealing rings for jointed glassware

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Joint grease, K.W.S

BUDDE - joint grease, silicone-free

Engler distilling flasks, DURAN®, ASTM, with beaded rim

Gas leak detection fluid, bottle