

Dropping funnels, cylindrical, with pressure equalizing tube, borosilicate glass 3.3, v

Cylindrical shape, with pressure equalizing tube. With NS socket and NS cone, plug with retaining device.

Material stopcock:

Glass key

With pressure compensation:

Yes

Nominal capacity ml	Division ml	Cone[s]	Socket[s]	Bore mm	Material	Shape	PK	Cat. No.
25	0.5	NS14/23	NS14/23	2.5	Boro 3.3	Cylindrical	1	6.072 172
50	1	NS14/23	NS14/23	2.5	Boro 3.3	Cylindrical	1	9.012 755
50	1	NS29/32	NS29/32	2.5	Boro 3.3	Cylindrical	1	9.012 767
100	2	NS14/23	NS14/23	2.5	Boro 3.3	Cylindrical	1	9.012 756
100	2	NS29/32	NS29/32	2.5	Boro 3.3	Cylindrical	1	9.012 768
250	5	NS14/23	NS14/23	4	Boro 3.3	Cylindrical	1	9.012 757
250	5	NS29/32	NS29/32	4	Boro 3.3	Cylindrical	1	9.012 769
500	10	NS29/32	NS29/32	4	Boro 3.3	Cylindrical	1	9.012 770
1000	20	NS29/32	NS29/32	6	Boro 3.3	Cylindrical	1	7.601 147

Zubehör:

PTFE sealing rings for jointed glassware

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Joint grease, K.W.S

BUDDE - joint grease, silicone-free

Engler distilling flasks, DURAN®, ASTM, with beaded rim

Gas leak detection fluid, bottle

Alternativen:

Dropping funnels, cylindrical, with pressure equalizing tube, borosilicate glass 3.3, with PTFE key