

Coil condensers, DURAN[®] tubing, with plastic olives

Made of DURAN[®] tubing. Acc. to DIN 12593 and in addition to DIN standard. With jacketed cooling coil. Water connection with GL 14 thread with two threaded plastic olives.

Connection size:

GL 14

Length jacket mm	Ground size	Cooling surface cm ²	Connection	Standard	Material	Type of construction	PK	Cat. No.
160	NS14/23	240	Tubing olive		DURAN [®] tubing	Intensive condenser	1	6.900 522
160	NS29/32	240	Tubing olive	DIN 12593	DURAN [®] tubing	Intensive condenser	1	9.012 542
250	NS24/29	400	Tubing olive		DURAN [®] tubing	Intensive condenser	1	6.310 602
250	NS29/32	400	Tubing olive	DIN 12593	DURAN [®] tubing	Intensive condenser	1	9.012 543
400	NS29/32	650	Tubing olive	DIN 12593	DURAN [®] tubing	Intensive condenser	1	9.012 544
500	NS29/32	840	Tubing olive		DURAN [®] tubing	Intensive condenser	1	6.203 018

Zubehör:

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Gas leak detection fluid, bottle

Joint grease, K.W.S

Engler distilling flasks, DURAN[®], ASTM, with beaded rim

BUDDE - joint grease, silicone-free

PTFE sealing rings for jointed glassware