

Spiral condensers, borosilicate glass 3.3, with PP olive

Made of borosilicate glass 3.3 which is resistant to heat and almost all chemicals. They meet the highest quality standards according to DIN 12592. The Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing. The side arms accept tubings with 8/9 mm bore size.

Standard:

DIN 12592

Length jacket mm	Ground size	Connection	Material	Type of construction	PK	Cat. No.
250	NS29/32	PP olive	Boro 3.3	Spiral condenser	1	6.257 820
400	NS29/32	PP olive	Boro 3.3	Spiral condenser	1	4.008 414

Zubehör:

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Gas leak detection fluid, bottle

Joint grease, K.W.S

Engler distilling flasks, DURAN®, ASTM, with beaded rim

BUDDE - joint grease, silicone-free

PTFE sealing rings for jointed glassware