

Receiver adapters, with olives, DURAN[®] tubing, bent 105°;

Made of DURAN[®] tubing. Acc. to DIN 12594, with NS socket and cone, with drain tip, vacuum connection: glass olive.

Connection: Glass olive
Standard: DIN 12594
Shape: Bent, 105°;

Cone[s]	Socket[s]	Material	PK	Cat. No.
NS14/23	NS14/23	DURAN [®] tubing	1	9.012 613
NS29/32	NS29/32	DURAN [®] tubing	1	9.012 623

Zubehör:

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Gas leak detection fluid, bottle

Joint grease, K.W.S

Engler distilling flasks, DURAN[®], ASTM, with beaded rim

BUDDE - joint grease, silicone-free

PTFE sealing rings for jointed glassware