

# Condensers, ground glass joint, Liebig, borosilicate glass 3.3, with PP olives

Made of borosilicate glass 3.3 which is resistant to heat and almost all chemicals. They meet the highest quality standards according to DIN 12576. Polypropylene side arm has a screw thread connector for easy and safe fitting of the tubing. The side arms accept tubings with 8/9mm bore size.

Standard:

DIN 12576

Length jacket mm	Ground size	Connection	Material	Type of construction	PK	Cat. No.
250	NS29/32	PP olive	Boro 3.3	Liebig condenser	1	9.012 546
400	NS29/32	PP olive	Boro 3.3	Liebig condenser	1	4.008 405

Zubehör:

PTFE sealing rings for jointed glassware

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Joint grease, K.W.S

BUDDE - joint grease, silicone-free

Engler distilling flasks, DURAN®, ASTM, with beaded rim

Gas leak detection fluid, bottle

Alternativen:

Condensers acc. to Allihn with GL threads, DURAN® tubing

Condensers acc. to Dimroth, with GL threads, DURAN® tubing