

Spiral condensers, DURAN[®] tubing, with PP olives

Made of DURAN[®] tubing. For cooling condensates. Water connections with GL14 threads with two threaded plastic olives.

Connection size: GL 14

Length jacket mm	Ground size	Cooling surface cm ²	Connection	Material	Type of construction	PK	Cat. No.
160	NS14/23	230	PP olive	DURAN [®] tubing	Spiral condenser	1	9.012 522
250	NS29/32	500	PP olive	DURAN [®] tubing	Spiral condenser	1	9.012 523
400	NS29/32	820	PP olive	DURAN [®] tubing	Spiral condenser	1	6.225 557
500	NS29/32	1040	PP olive	DURAN [®] tubing	Spiral condenser	1	9.012 524

Zubehör:

PTFE sealing rings for jointed glassware

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Joint grease, K.W.S

BUDDE - joint grease, silicone-free

Engler distilling flasks, DURAN[®], ASTM, with beaded rim

Gas leak detection fluid, bottle