

Drying tubes, U-shaped, DURAN^{®} tubing, with interchangeable ground joint stopcock

Made of DURAN^{®} tubing. DIN 12616. U-shape, side tubes plain or with olives (outer dia. 6 mm).

Shape: U-shape
With stopcock: Yes
Connection: Tube connection
Standard: DIN 12616

Length mm	Socket[s]	Material	PK	Cat. No.
125	NS14/23	DURAN ^{&reg;} tubing	1	6.902 239
150	NS19/26	DURAN ^{&reg;} tubing	1	9.253 115

Zubehör:

Sleeves, cone joint, PTFE
Sleeves with ribs, PTFE
Joint clips for conical joint sleeve connections, stainless steel, with set screw
Ground joint clips for sleeve connections, wire, Chrome-nickel steel
Joint clips KECK, POM, for conical ground joints
Joint clips KECK, POM, for spherical ground joints
Gas leak detection fluid, bottle
Joint grease, K.W.S
Engler distilling flasks, DURAN[®], ASTM, with beaded rim
BUDDE - joint grease, silicone-free
PTFE sealing rings for jointed glassware

Alternativen:

Drying tubes, DURAN[®] tubing, straight
Drying tubes, PP