

Kjeldahl flasks, DURAN^{®}

- DIN 12360
- With beaded rim

Nominal capacity ml	Neck Ø mm	Ø mm	Height mm	With beaded rim	Material	PK	Cat. No.
50	22	51	200	Yes	Boro 3.3	1	9.141 817
100	22	60	200	Yes	Boro 3.3	1	9.141 824
250	34	81	270	Yes	Boro 3.3	1	9.141 836
500	34	101	300	Yes	Boro 3.3	1	9.141 844
750	34	115	340	Yes	Boro 3.3	1	9.141 851
1000	34	126	350	Yes	Boro 3.3	1	9.141 854

Note: special sizes available on request.

Zubehör:
Flask rings, pressed cork
Flask support rings, "BiBase", silicone elastomer
Round bottom flask holder, PP