

# Rectangular retort stand bases, enamel, coated steel plate

Highly stable laboratory stand due to the steel rectangular base which is stove enamel coated. With an M10 threaded hole in the middle of one short side. Rubber support feet prevents the stand from sliding. Laboratory stands are supplied without support rod. Please select the suitable rod from "support rods".

Material: Steel  
Material coating: Enamel

Width mm	Depth mm	Dimensions (WxD)	Weight g	PK	Cat. No.
130	210	130 x 210 mm	1725	1	6.242 701
150	300	150 x 300 mm	2830	1	4.008 486
200	315	200 x 315 mm	4000	1	4.008 487

Zubehör:  
Support rings  
Support ring  
Support rods, stainless steel  
Support rods galvaniser steel

Alternativen:  
Retort stand bases  
Retort stand bases