

# Condensers acc. to Liebig, DURAN<sup>&reg;</sup> tubing, with plastic olives

Made of DURAN<sup>&reg;</sup> tubing. DIN 12576 and in addition to DIN standard.

With NS cone and socket. Water connection GL14 threads with two threaded plastic olives.

Connection size:

GL 14

Standard:

DIN 12576

Length jacket mm	Ground size	Cooling surface cm <sup>2</sup>	Connection	Material	Type of construction	PK	Cat. No.
160	NS14/23	60	Tube connection	DURAN <sup>&amp;reg;</sup> tubing	Liebig condenser	1	9.012 502
160	NS19/26	90	Tube connection	DURAN <sup>&amp;reg;</sup> tubing	Liebig condenser	1	9.012 501
250	NS14/23	100	Tube connection	DURAN <sup>&amp;reg;</sup> tubing	Liebig condenser	1	9.012 503
250	NS29/32	130	Tube connection	DURAN <sup>&amp;reg;</sup> tubing	Liebig condenser	1	9.012 505
400	NS29/32	210	Tube connection	DURAN <sup>&amp;reg;</sup> tubing	Liebig condenser	1	9.012 506

## Zubehör:

PTFE sealing rings for jointed glassware

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Joint grease, K.W.S

BUDDE - joint grease, silicone-free

Engler distilling flasks, DURAN<sup>®</sup>, ASTM, with beaded rim

Gas leak detection fluid, bottle

## Alternativen:

Condensers acc. to Allihn with GL threads, DURAN<sup>®</sup> tubing

Condensers acc. to Dimroth, with GL threads, DURAN<sup>®</sup> tubing