

Dropping funnels, cylindrical, without pressure equalizing tube, Borosilicate glass 3

PTFE stopcock, PE stoppers. Cylindrical shape, without pressure equalising tube. ISO 4800 and DIN 12242.

Material stopcock:

PTFE key

Standard:

ISO 4800, DIN 12242

With pressure compensation:

No

Nominal capacity ml	Cone[s]	Socket[s]	Material	Shape	PK	Cat. No.
100	NS29/32	NS19/26	Boro 3.3	Cylindrical	1	9.140 411
250	NS29/32	NS29/32	Boro 3.3	Cylindrical	1	9.140 412
500	NS29/32	NS29/32	Boro 3.3	Cylindrical	1	9.140 413
1000	NS29/32	NS29/32	Boro 3.3	Cylindrical	1	9.140 414

Zubehör:

Sleeves, cone joint, PTFE

Sleeves with ribs, PTFE

Joint clips for conical joint sleeve connections, stainless steel, with set screw

Ground joint clips for sleeve connections, wire, Chrome-nickel steel

Joint clips KECK, POM, for conical ground joints

Joint clips KECK, POM, for spherical ground joints

Gas leak detection fluid, bottle

Joint grease, K.W.S

Engler distilling flasks, DURAN®, ASTM, with beaded rim

BUDDE - joint grease, silicone-free

PTFE sealing rings for jointed glassware