

Water baths PURA™

Flexible in use for all purposes in research and quality control. All PURA™ water baths (except PURA™ 4) have an integrated drain screw. The bright, easy-to-read white LED display is clearly visible, even from a distance. With integrated side recessed grips for easy transport and non-slip rubber feet for safe standing. The complete interior surface of the water baths can be used, as there are no installation elements in the interior.

- Including platform
- Integrated timer
- Bright display with membrane keypad
- Splash-proof
- Built-in dry-running protection

Please order bath cover separately.

Temperature: 18 ... 99,9 °C
 Temperature stability: ±0,15 °C
 Plug type: EU
 Power supply: 230 V, 50/60 Hz
 IP code: IP 21

Type	Min. capacity l	Max. capacity l	Capacity	Min. temperature °C	Max. temperature °C	Display	Type of display	With ci
PURA™ 4	0.8	4.5	0.8 ... 4.5 l	18	99,9	Digital	LED	No
PURA™ 10	1.4	9.5	1.4 ... 9.5 l	18	99,9	Digital	LED	No
PURA™ 14	2	14	2.0 ... 14.0 l	18	99,9	Digital	LED	No
PURA™ 22	3.4	25.5	3.4 ... 25.5 l	18	99,9	Digital	LED	No
PURA™ 30	4.8	36	4.8 ... 36.0 l	18	99,9	Digital	LED	No

Zubehör:

- Hinged bath covers for water baths PURA™
- Cooling coil for water baths PURA™
- Accessories for Water baths, PURA™
- Test tube racks for water baths PURA™, shaking water baths SW and heating circulators CORIO™
- Liquid level/cooling set for water baths PURA™ and shaking water baths SW
- Flat bath covers with ring sets for water baths PURA™